

## Mating Protocol

Preliminary calculations:

Want 250 $\mu$ l of OD600=10 per filter, or 2.5 OD's/filter for EACH *E.coli* strain.

Want 500 $\mu$ l of 100 $\mu$ g/ml chla for each filter, or 50 $\mu$ g Chla/filter.

1. Have 500ml exponential ammonia grown *Nostoc* ready. Take chla reading and harvest enough for all matings. Resuspend to 100 $\mu$ g chla/ml and use 500 $\mu$ l per filter.
2. Inoculate 1/40th overnight culture of each *E.coli* strain into selective LB. Check OD600 and use when near 1. Calculate total OD600 units per culture. Centrifuge and resuspend to OD=10. Mix *E. coli* strains 1:1, to allow for 500 $\mu$ l per mating filter. Centrifuge and wash twice in LB.
3. Wash Millipore HATF nitrocellulose filters in d water and place between Watman #1 filters. Autoclave in glass petri dish dry cycle( 20-20min). Place atop AA mops, ammonia, 0.5% LB mating plates.
4. Mix 500 $\mu$ l each *Nostoc* and *E.coli* from step 1&2 in a microfuge tube and spin briefly in microfuge. Pipet off all but ~200-300 $\mu$ l of supernatant and resuspend by pipetting up and down. Place atop mating filter (step 3) and spread. Allow for moisture to be absorbed/dried before parafilming.
5. Place in low light (bottom shelf CO2 incubator) overnight. Transfer filter to AA mops ammonia and allow to recover in low light. When cells look healthy transfer filter to AA, mops ammonia, antibiotic plates, and place in high light